

# STREAMLINED POSTERIOR RESTORATION PROCEDURE USING UNIVERSAL COMPOSITE

with CLEARFIL MAJESTY™ ES-2 Universal



CASE BY  
DANIEL VASQUEZ, DDS

The dental practitioner's time is the most valuable factor in the dental practice. Hence, it should be well spent, and saved wherever possible. In the context of placing direct posterior restorations, the universal composite CLEARFIL MAJESTY™ ES-2 Universal provides valuable support to anyone who would like to achieve this goal - being universally applicable, it allows users to spend less time on material selection. With its universal shade concept consisting of a single posterior shade and two anterior shade options, it eliminates the need for shade determination. This gives users more time to focus on fulfilling their patient's needs. The following clinical case illustrates the use of the innovative material in the posterior region.



Fig. 1.

Pre-operative situation with multiple carious lesions in the second premolar and both molars.



Fig. 2.


Isolated working field.



**Fig. 3.**  
View of the quadrant after cavity preparation.



**Fig. 4.**  
Selective etching of the enamel with phosphoric acid gel.



**Fig. 5.**  
Dispensing of the universal adhesive CLEARFIL™ Universal Bond Quick into a mixing dish.



**Fig. 6.**  
Application of the universal adhesive to the enamel and dentin after etching, rinsing and air-drying.



**Fig. 7.**  
Lining up of the cavity floors with of a thin layer of flowable composite (CLEARFIL MAJESTY™ Flow).



**Fig. 8.**  
Immediate treatment outcome after the placement of several layers of CLEARFIL MAJESTY™ ES-2 Universal in the shade Universal.



## FINAL SITUATION

**Fig. 9.**  
Post-operative image taken after rubber dam removal.

### **CLEARFIL MAJESTY™ ES-2 Universal**

VISIT THE PRODUCT PAGE ONLINE FOR MORE DETAILS

#### **WHAT CAN YOU EXPECT TO FIND ONLINE**

- ✓ General information about CLEARFIL MAJESTY™ ES-2 Universal
- ✓ Full Product Assortment
- ✓ Safety Data Sheets



**KURARAYNORITAKE.EU**