

Which Sonicare interdental cleaner is right for your patient?



Power Flosser

For patients who:

- 1 Practice good oral care, but struggle with flossing
- 2 Would benefit from a more thorough, deep clean
- 3 Are looking for a more customizable power flossing experience




Cordless Power Flosser

For patients who:

- 1 Are starting a regular flossing habit
- 2 Seek convenience and value portability
- 3 Have more limited counter space



A thorough clean, effortlessly



Quad Stream

Unique cross-shaped nozzle creates four wide streams for a thorough clean interproximally and along the gingival margin.



Pulse Wave technology

Brief pauses in the cleaning pulsations prompt your patients to glide the nozzle to the next interproximal space to help ensure a whole-mouth clean.

Call (800) 422-9448 or visit philips.com/cordlesspowerflosser to learn more



© 2022 Koninklijke Philips N.V. (KPNV). All rights reserved. PHILIPS and the Philips shield are trademarks of KPNV. SONICARE and the Sonicare logo are trademarks of KPNV. Other trademarks and trade names are those of their respective owners.

www.philips.com
7510869

PHILIPS

sonicare

Interdental cleaning



Flossing reinvented

Introducing the Philips Sonicare Power Flossers

Advanced, effective interproximal cleaning solutions

Philips Sonicare Cordless Power Flosser

The easy way to clean in between



Sleek ergonomic design and 360-degree rotating nozzle allow for easy maneuverability



Three intensity settings and two flossing modes: Clean and Deep Clean¹



Removes up to 99.9% of biofilm in treated areas²



Up to 150% more effective than floss in improving gum health³



Philips Sonicare Power Flosser

A thorough clean, effortlessly



10 intensity settings and two flossing modes for a highly personalized flossing experience



Whisper quiet technology allows for quiet operation without compromising effectiveness



Up to 180% more effective than string floss at improving gingival health⁵

1 Features vary by model
2 In an in-vitro study, actual in-mouth results may vary
3 Between teeth, after two weeks, in patients with moderate to severe gingivitis, when used with the Quad Stream nozzle on high setting with a manual toothbrush vs. using a manual toothbrush with string floss
4 Vs leading competitor, in an in-vitro study, results may vary
5 When used with the Quad Stream nozzle on setting 8 with a manual toothbrush in patients with moderate to severe gingivitis vs. using a manual toothbrush with string floss