

# 10 Keys to Perfect IOS Cases

Master Digital Impressions & Elevate Your Lab Results



## 3 Biggest Mistakes that Lead to Remakes

### **QUAD SCANS INSTEAD OF FULL ARCH**

Partial scans increase likelihood of high occlusion and open bites, as well as inability to match contralateral contours.

### **UNCLEAR OR INCOMPLETE MARGINS**

The #1 reason for remakes – if we can't see the margin on our screen, our crown margin is an educated guess at best.

### **VAGUE PRESCRIPTIONS**

Missing patient information could lead to guesswork, not precision.

## 7 Actionable Tips for Flawless Digital Impressions

- 1. USE DOUBLE-PACK RETRACTION CORD**  
Creates ideal margin visibility – the foundation of perfect restorations.
- 2. ENSURE SCANS ARE DEBRIS-FREE**  
Digital scanners capture EVERYTHING – including blood, moisture, and debris.
- 3. REVIEW YOUR SCANS BEFORE SUBMISSION**  
Take an extra minute to check for missing areas, distortions, or artifacts.
- 4. SUBMIT FULL ARCH SCANS WHEN POSSIBLE**  
Provides complete context for proper fit, alignment, and function.
- 5. AVOID MULTIPLE SCANS OF THE SAME AREA**  
Creates overlapping data that complicates processing.
- 6. SUBMIT SCANS IMMEDIATELY**  
Keeps workflow moving and prevents delays.
- 7. PROVIDE COMPREHENSIVE PATIENT DETAILS**  
Note grinding habits, aesthetic preferences, and specific challenges.

## The #1 Success Factor: Communication

The most successful cases come from doctors who provide more information than necessary, anticipate potential challenges, and maintain open channels for collaboration.

### **CONTACT US**

Partner with our master ceramists and certified dental technicians for exceptional results

