



More Patients Than Ever Can Now Be Treated with LOCATOR®

INTRODUCING LOCATOR ANGLED ABUTMENT

Synergy Members
Receive 20% Off
When Ordering At
Zestdent.com

The new LOCATOR Angled Abutment connects clinicians to a greater number of edentulous patients by offering an affordable, simple and flexible solution that meets their unique dental healthcare needs. This innovation adds an additional 15 degrees of angle correction and makes it possible to treat more patients than ever in combination with the LOCATOR Removable or the LOCATOR FIXED restorative systems.

> DIVERGENCE CORRECTION

LOCATOR Angled Abutments correct implant angulation by 15 degrees. The LOCATOR system continues to correct up to an additional 20 degrees of divergence with the use of LOCATOR removable extended range or LOCATOR FIXED inserts.

> REMOVABLE OR FIXED

Compatible with LOCATOR Removable Extended Range inserts and LOCATOR FIXED Inserts.

> TREATMENT FLEXIBILITY

Angled abutments offer greater opportunity to avoid critical anatomy in implant placement and increase A-P spread to limit distal cantilevers when using LOCATOR FIXED.

> SIMPLE SOLUTION THAT REDUCES COMPLICATIONS

Properly aligned abutments not only lower risks of fractures and loosening but also ensures better retention and long-term stability. Use of this standard abutment also avoids the need for custom abutments.

> SELF-ALIGNING DESIGN

The self-aligning feature of the LOCATOR attachment allows a patient to easily seat their overdenture without the need for accurate alignment, and without causing damage to the attachment components. This self-aligning feature also increases the longevity of the LOCATOR abutment.

> TiN COATING

Features the same LOCATOR TiN coating you know and trust! LOCATOR Angled abutments allow a better path of insertion in highly angulated cases, even further protecting the abutment from wear.



REMOVABLE

*Same
LOCATOR
Abutment,
Two Solutions*



FIXED



LOCATOR Angled Abutment Ordering Information

CONNECTION	PART NUMBER	CUFF HEIGHT	REPLACEMENT SCREW	REMOVAL TOOL
Nobel NP	02415-SB	2.5mm	02414-SB	02412-SB
	02416-SB	3.5mm		
	02417-SB	4.5mm		
	02418-SB	5.5mm		
	02419-SB	6.5mm		
	02420-SB	7.5mm		
Nobel RP	02422-SB	2.5mm	02421-SB	02413-SB
	02423-SB	3.5mm		
	02424-SB	4.5mm		
	02425-SB	5.5mm		
	02426-SB	6.5mm		
	02427-SB	7.5mm		
ZimVie TSV 3.5	02429-SB	2.5mm	02428-SB	N/A
	02430-SB	3.5mm		
	02431-SB	4.5mm		
	02432-SB	5.5mm		
	02433-SB	6.5mm		
	02434-SB	7.5mm		
ZimVie TSV 4.5	02435-SB	2.5mm	02428-SB	N/A
	02436-SB	3.5mm		
	02437-SB	4.5mm		
	02438-SB	5.5mm		
	02439-SB	6.5mm		
	02440-SB	7.5mm		
Neodent GM	02442-SB	2.5mm	02441-SB	N/A
	02443-SB	3.5mm		
	02444-SB	4.5mm		
	02445-SB	5.5mm		
	02446-SB	6.5mm		
	02447-SB	7.5mm		
Straumann BLX	02398A-SB	2.5mm	02404-SB	02412-SB
	02399A-SB	3.5mm		
	02400A-SB	4.5mm		
	02401A-SB	5.5mm		
	02402A-SB	6.5mm		
	02403A-SB	7.5mm		

Required Tooling

TOOL	PART NUMBER
LOCATOR Angled Abutment Driver Set	02405-SB
LOCATOR Enhanced Core Tool	08340-M
LOCATOR Angled Abutment Removal Tool*	02413-SB or 02412-SB

*Check connection above for Removal Tool requirements