

Icon Smooth Surface



The only microinvasively **permanent esthetic solution** to white spots

Icon is the only microinvasive treatment solution proven to diminish the appearance of white spot enamel lesions, up to D1, in one patient visit.

Types of treatable lesions include:



Post orthodontic white spots



Mild-to-moderate fluorosis



Congenital enamel defects



- ✓ **No** anesthesia
- ✓ **No** drilling
- ✓ **One** visit

Developmental lesions

FEATURES AND BENEFITS

- Esthetic results on smooth surfaces; diminishes the appearance of white spots
- Caries are arrested at an early stage
- Healthy tooth structure is preserved
- No anesthesia or drilling required
- Teeth are restored back to their natural shade
- One quick and easy patient visit
- Available insurance code D2990



The proven microinvasive solution for white spots and proximal caries

Revolutionary caries infiltration with Icon fills the gap between prophylaxis and filling. It also optically blends in the lesion with the healthy enamel, for esthetic treatments of front teeth, where appearance plays an important role.

Ingeniously simple

The Icon infiltration method is simple. After pretreatment with Icon Etch, a highly fluid resin (Icon Infiltrant) is applied to the affected area. Capillary action draws the infiltrant deep into the porous enamel, and it is then light-cured.

Immediately following treatment, the infiltrated area visually takes on the appearance of the surrounding healthy enamel, and the patient leaves with a transformative esthetic result.

Healthy tooth structure is preserved.



The light refraction of the otherwise translucent enamel is changed and the lesion has a whitish-opaque appearance.

The infiltration fills the porous tissue while the lesion takes on the appearance of the surrounding healthy enamel.

220402

Icon Smooth Surface Mini Kit

2 Patient Packs each including:
1-.45ml Syringe Icon Etch
1-.45ml Syringe Icon Dry
1-.45ml Syringe Icon Infiltrant
6 Smooth Surface Tips
1 Luer-Lock Tip

220403

Icon Smooth Surface Cube

7 Patient Packs each including:
1-.45ml Syringe Icon Etch
1-.45ml Syringe Icon Dry
1-.45ml Syringe Icon Infiltrant
6 Smooth Surface Tips
1 Luer-Lock Tip

220385

Icon Etch Refill

3-.45ml Syringes Icon Etch
15 Smooth Surface Tips

ICON AND WHITE SPOTS



Initial situation: White spots after removing braces.



Result after treating with Icon Etch.



Drying and visual check with Icon Dry.



Application of Icon Infiltrant.



Repeated application of Infiltrant.



Result directly after light-curing.



Result after polishing.



One week after treatment.



Scan to learn more about Icon Smooth Surface and Icon Proximal, order waiting room displays, and become an Icon Provider!

Plus, learn about Icon Proximal, proven to effectively treat interproximal caries up to D1.

