


INDIRECT RESTORATIVE MATERIALS

ZIRCONIA - Monolithic or Layered

Bruxzir

Crown Removal



Endo Access



Removal	Grit	
Great White®Z 856-018M	M	
Great White®Z 850-018M	M	
Great White®Z 801-016M	M	
Great White®Z 801-018M	M	

Preparation	Grit	
Piranha® 856-024C (mesial buccal distal)	C	
Piranha® 379-023SC (occlusal surface)	SC	
Revelation® 878K-018SC	SC	
Revelation® 811-037SC	SC	
Great White® Ultra Carbide 856-018		
Great White® Ultra Carbide 856-016		
Great White® Ultra Carbide 856-020		
Great White® Ultra Carbide 379-023 (occlusal surface)		

Adjust/Contour	Grit	
Revelation® 368-020F (occlusal)	F	
Revelation® 856-021F (proximal)	F	
Revelation® 858-014F (marginal)	F	

Polish	Grit	
Jazz® polisher Medium Cup ZA2S	M	
Jazz® polisher Fine Cup ZA2S	F	
Jazz® polisher Medium Flame ZA2S	M	
Jazz® polisher Fine Flame ZA2S	F	

PORCELAIN/MILLED CERAMICS

Feldspathic leucite reinforced porcelain (Empress) - Lithium disilicate porcelain crowns (Emax) - LAVA or Procera

Crown Removal



Endo Access

Removal	Grit	
Piranha® 856-016C	M	
Piranha® 856-016C	M	
Great White®Z 801-016M	M	
Great White®Z 801-018M	M	

Preparation	Grit	
Piranha® 856-024C (mesial buccal distal)	C	
Piranha® 379-023SC occlusal surface	SC	
Revelation® 878K-018SC	SC	
Revelation® 811-037SC	SC	
Great White® Ultra Carbide 856-018		
Great White® Ultra Carbide 856-016		
Great White® Ultra Carbide 856-020		
Great White® Ultra Carbide 379-023 (occlusal surface)		

Adjust/Contour	Grit	
Revelation® 368-020F (occlusal)	F	
Revelation® 856-021F (proximal)	F	
Revelation® 858-014F-S (marginal)	F	
Jazz® polisher Coarse Cup P3S	C	
Jazz® polisher Coarse Cup P3S	C	

Polish	Grit	
Jazz® polisher Medium Cup P3S	M	
Jazz® polisher Fine Cup P3S	F	
Jazz® polisher Medium Flame P3S	M	
Jazz® polisher Fine Flame P3S	F	

PORCELAIN (Fused-to Metal)

Ceramic Portion

Metal Portion

Removal	Grit	
Piranha® 856-016C	M	
Piranha® 856-016C	M	
FG Carbide 1558		
FG Carbide 332L		
FG Carbide 558S		
FG Carbide702		
GW2 (1558)		

Preparation	Grit	
Piranha® 856-024C (mesial buccal distal)	C	
Piranha® 379-023SC occlusal surface	SC	
Revelation® 878K-018SC	SC	
Revelation® 811-037SC	SC	
Great White® Ultra Carbide 856-018		
Great White® Ultra Carbide 856-016		
Great White® Ultra Carbide 856-020		
Great White® Ultra Carbide 379-023 (occlusal surface)		

Adjust/Contour	Grit	
Revelation® 368-020F (occlusal)	F	
Revelation® 856-021F (proximal)	F	
Revelation® 858-014F-S (marginal)	F	
Jazz® polisher Coarse Cup P3S	C	
Jazz® polisher Coarse Cup P3S	C	

Polish	Grit	
Jazz® polisher Medium Cup P3S	M	
Jazz® polisher Fine Cup P3S	F	
Jazz® polisher Medium Flame P3S	M	
Jazz® polisher Fine Flame P3S	F	

GOLD ALLOYS (Nobel)

FG Carbide 1558		
FG Carbide 332L		
FG Carbide 558		
FG Carbide 702		

Piranha® 856-024C-SS	C	
Piranha® 379-023SC	SC	
Revelation® 878K-018SC	SC	
Revelation® 811-037SC	SC	
Great White® Ultra Carbide 856-018		
Great White® Ultra Carbide 856-016		
Great White® Ultra Carbide 856-020		
Great White® Ultra Carbide 379-023		

Revelation® 368-020F (occlusal)	F	
Revelation® 856-021F (proximal)	F	
Revelation® 858-014F (marginal)	F	
Jazz® polisher Coarse Cup P3S	C	
Jazz® polisher Coarse Cup P3S	C	

Jazz® polisher Medium Cup P3S	M	
Jazz® polisher Fine Cup P3S	F	
Jazz® polisher Medium Flame P3S	M	
Jazz® polisher Fine Flame P3S	F	

NICKEL OR COBALT-CHROME (Base-Metal) Alloys or Titanium Abutments

FG Carbide 1558		
FG Carbide 332L		
FG Carbide 558		
FG Carbide 702		

Piranha® 856-024C	C	
Piranha® 379-023SC	SC	
Revelation® 878K-018SC	SC	
Revelation® 811-037SC	SC	
Great White® Ultra Carbide 856-018		
Great White® Ultra Carbide 856-016		
Great White® Ultra Carbide 856-020		
Great White® Ultra Carbide 379-023		

Revelation® 368-020F (occlusal)	F	
Revelation® 856-021F (proximal)	F	
Revelation® 858-014F (marginal)	F	
Jazz® polisher Coarse Cup P3S	C	
Jazz® polisher Coarse Cup P3S	C	

Jazz® polisher Medium Cup P3S	M	
Jazz® polisher Fine Cup P3S	F	
Jazz® polisher Medium Flame P3S	M	
Jazz® polisher Fine Flame P3S	F	

DIRECT RESTORATIVE MATERIALS

AMALGAM

Removal		Preparation		Adjust/Contour	Grit		Polish	Grit	
Carbide 1558		Carbide 1558		Revelation® 368-020F (occlusal)	F		Jazz® polisher Medium Cup P3S	M	
Carbide 332L		Carbide 332L		Revelation® 856-021F (proximal)	F		Jazz® polisher Fine Cup P3S	F	
Carbide 558S		Carbide 558		Revelation®858-014F (marginal)	F		Jazz® polisher Medium Flame P3S	M	
Carbide 702		Carbide 702		Jazz® polisher Coarse Cup P3S	C		Jazz® polisher Fine Flame P3S	F	
				Jazz® polisher Coarse Cup P3S	C				

DIRECT AND INDIRECT RESTORATIVE MATERIALS

- COMPOSITE
- GLASS IONOMER restoration
- RESIN IONOMER restoration

Removal		Preparation		Adjust/Contour	Grit		Polish	Grit	
Carbide 1558		Carbide 1558		Revelation® 368-020F (occlusal)	F		Jazz® polisher Fine Cup C1S	F	
Carbide 332L		Carbide 332L		Revelation® 856-021F (proximal) & Revelation® 858-014F (marginal)	F		Jazz® polisher Medium Flame C1S	F	
Carbide 558		Carbide 558							
Carbide 702		Carbide 702							

Better Patient Outcomes

Improved Efficiency

Faster Practice Growth



1145 Towbin Avenue • Lakewood, NJ
sswhitedental.com

#64669

Dentists Materials Reference