



Accur8[™] & Matrix1[™] IMPRESSION MATERIALS

ACCUR8™

STARTER BUNDLE

- (2) 024425-BAccur8 Light Body Fast Set* 2PK
- (2) 022405-B Accur8 Heavy Body Fast Set* 2PK
- (2) 023405-B Accur8 Rigid Fast Set* 2PK
- (1) 026400-B Bite Registration 2PK
- (1) 047501-B Accur8 Dual-Arch Impression Trays 50PK

Retail Value: \$649.64 Bundle Price: \$341.02 *Mix & Match*

MATRIX1[™]

- (2) 018902-B Matrix1 Material 2PK
- (1) K0181 Spear Provisional System
- (1) 046050-B Posterior Impression Tray 50PK
- (1) 046150-B Anterior Impression Tray 50PK

Retail Value: \$591.38 Bundle Price: \$290.81

Make a Great Impression with Your Patients

With Accur8[™] and Matrix1[™] Impression Materials you can rely on product solutions that produce the results you expect—consistently. *Substitute the Regular Set for Fast Set

Impression materials for every procedure and case

Impression Material		Packaging	Order No.	Intended Use	Working Time (min)	Min. Set Time (min)	Max. Total Cure Time (min)
Accur8 Light Body High Flow	Regular Set	2-Pack	024420-В	Flows on its own under the hydraulic pressure of the tray material and flows effortlessly into the sulcus.	2:00	3:30	5:30
		12-Pack	019522-B		2.00	5.50	0.00
	Fast Set	2-Pack	024425-B		1:00	1:45	2:45
		12-Pack	019525-B				
Accur8 Heavy Body High Flow	Regular Set	2-Pack	023400-В	For use on full-arch impressions. Blends with the light body, preventing over displacement. Cures to an elastic 68-70 durometer. Can be used for aligner impressions.	2:00	3:30	5:30
		12-Pack	013912-В				
	Fast Set	2-Pack	023405-B		1:00	1:45	2:45
		12-Pack	013952-B				
	Regular	2-Pack	022400-B			5:30	
Accur8 Rigid Tray High Flow	Set	12-Pack	012912-B	Specifically designed for use with a dual-arch impression technique. A unique low viscosity combined with high rigidity when set. Cures to a very rigid 80 durometer. Provides more rigidity than a putty or a heavy body. Designed to prevent impression distortion.		5.50	5.50
	Fast Set	2-Pack	022405-B		1:00	1:45	2:45
		12-Pack	012952-B				
Accur8 Bite Registration		2-Pack	026400-В	Mousse-like bite registration, not brittle, with an ideal flow to capture accurate detail. Carvable, ultra-high 90 durometer cure eliminates distortion during prosthetic fabrication. Records relationship between the upper and lower teeth when jaw is in a closed position.	0:30	0:45	1:15
		12-Pack	016912-В				
Matrix1		2-Pack	018902-B	Ultra-quick, ultra-accurate silicone matrix material used to fabricate a matrix for provisionals and to form a lingual matrix.	0:30	0:30	1:00
WIDUIX1		12-Pack	018912-B				

35% OFF ANY QUANTITY

Accessories

Product Description	Order No.	Qty/Pk
Dual-Port 1:1 Impression Material Dispensing Gun	5027164U0	1
Accur8™ Dual-Arch Posterior Impression Trays	047501-B	50
Anterior Temporary Impression Trays for use with Matrix1™	046150-В	50
Posterior Temporary Impression Trays for use with Matrix1™	046050-В	50
6.5 Dual-Port Mixing Tips (Green) for use with Accur8™ Heavy Body, Rigid Tray, Bite Registration and Matrix1™	033526-В	50
6.5 Dual-Port Mixing Tips (Green) for use with Accur8™ Heavy Body, Rigid Tray, Bite Registration and Matrix1™	033126-В	100
4.2 Dual-Port Mixing Tips (Yellow) for use with Accur8™ Light Body	033521-B	50
4.2 Dual-Port Mixing Tips (Yellow) for use with Accur8™ Light Body	033121-B	100
Intraoral Tips (Yellow) for use with Accur8™ Light Body	032520-В	50
Intraoral Tips (Yellow) for use with Accur8™ Light Body	032120-В	100

Visit our website at BrasselerUSA.com To order call 800.841.4522 or fax 888.610.1937. ©2023 Brasseler USA.All rights reserved. Printed in the U.S.A.

f ♥ in ⊙ B-5759-05.23-5M-PRO