

# QUATTROcare™ Plus

Less time maintaining,  
more instrument sustaining

The user friendly QUATTROcare™ Plus handpiece maintenance system is designed to lubricate and purge your instruments with the simple push of a button.

Simple operation and loading make the QUATTROcare Plus easy to use. Maintenance with the QUATTROcare Plus system means less wear, less repairs, less staff down time, smoother operation, and increased longevity of your handpieces.

The QUATTROcare Plus can be used for virtually all manufacturer's high-speeds, low-speed contra-angles, straight handpieces and air driven scalers.



### Easy to use

Up to four (4) handpieces can be lubricated and purged in the same cycle.



### Safe, controlled purging

Completely controlled purging procedure.



### Economical investment

Can be used for virtually all manufacturers' handpieces, air driven scalers and air motors.



### Fast lubrication cycle

Just one minute per care cycle for four instruments, all programmatically and automatically.



### Chuck care maintenance

Chuck care system provides easy maintenance of handpiece chuck.



**QUATTROcare Plus**  
Part no. 1008.3805

# QUATTROcare™ Plus Maintenance System

## QUATTROcare Plus

Automatic handpiece maintenance and purging system. Correct maintenance of dental handpieces is critical to ensure top performance and longevity. KaVo QUATTROcare Plus is equipped with one MULTIflex adapter and three INTRAMatic (E-Type ISO 3964) adapters. Unit also includes one can of QUATTROcare Plus spray.

### Dimensions

- Height: 380 mm (1.25 feet), Width: 350 mm (1.15 feet), Depth: 230 mm (0.75 feet)
- Weight: 10 kg (22 pounds)
- Width with door ajar (right side - to exchange spray can): 16<sup>15</sup>/<sub>16</sub>". Depth with flap open (front): 18<sup>7</sup>/<sub>8</sub>"

### Technical

- Supply voltage 100 to 240 VAC
- Supply frequency 50 or 60 Hz
- Power consumption 41 VA
- Air consumption: 50 NL / min or 1.765 cfm
- Simply connects to ¼" airline
- 58 to 87 psi air pressure recommended (max 145 psi)
- Built in automatic pressure regulator controls flow

## Easy As...



## QUATTROcare Spray

In general, automated care of rotating dental handpieces using the KaVo QUATTROcare Plus is preferred over manual care. The handpieces are lubricated and serviced professionally and correctly with the proper amount of QUATTROcare Spray.

QUATTROcare Spray connects quickly and easily into the QUATTROcare Plus handpiece maintenance system.



### QUATTROcare Spray

Part no. 1005.4524