






CUT, FINISH & POLISH GUIDE

DIRECT & INDIRECT RESTORATIONS



BRASSELER
CANADA®

TABLE OF CONTENTS

Carbides

Amalgam Removal	4
Decay Removal	4
Preparation	4
Trimming and Finishing	4,8
Metal Cutting	4
Trimming Alloy, Acrylic, and Soft Silicone	11
Orthodontic Debonding	12
Surgical	12
Access	12

Diamonds

Amalgam Removal	4
Composite Removal	4
Preparation	4,8
Finishing	4,8
Porcelain Removal	6
Zirconia and Lithium Disilicate Removal	6
Depth Cutting	6,8
Breaking Contact	6
Reduction	6,7,8
Margin Refinement	7
Adjusting	9
Separating and Contouring Acrylic and Ceramics	11
Periodontal Ligament Separation	12
Bone Cutting	12
Pilot Access	12
Access	12
Sterile Preparation	13

Strips

Interproximal	5
Separating	5

Polishers & Abrasives

Single-Patient Composite	5,9
Multi-Use Composite	5,11
Final Polishing Composite	5
Shaping and Polishing Embrasures	5
Adjusting Ceramic Crowns	9
Polishing Zirconia	9,10
Polishing Lithium Disilicate	9,10
Polishing Feldspathic Porcelain	9,10
Adjusting All-Ceramics	10
Adjusting and Polishing Acrylic and Soft Silicone	11
Adjusting and Polishing PMMA	11

DIRECT RESTORATION

Brasseler's recommended diamonds, carbides, and polishers for each step of a Direct Restoration.

COMPOSITE

Amalgam Removal Diamond



Us No.	—
Order No.	6813
Head Size (1/10 mm)	014
Cutting Length (mm)	1.5
Grit	
Shank Type	FG

Amalgam Removal Carbides



Sabrecut®		
Us No.	—	1558SC
Order No.	H34	H031RS
Head Size (1/10 mm)	010	012
Cutting Length (mm)	2.0	4.0
Blade	—	—
Shank Type	FG	FG

Composite Removal Diamonds



ET® Series		
Us No.	—	DET6F
Order No.	6830	134F
Head Size (1/10 mm)	014	014
Cutting Length (mm)	2.7	6.0
Grit		
Shank Type	FG	FG

Decay Removal Carbide



Us No.	4
Order No.	H1
Head Size (1/10 mm)	014
Cutting Length (mm)	—
Blade	—
Shank Type	RA

Preparation Diamonds



Micro			
Us No.	330D	557D	—
Order No.	D232	D835	838M
Head Size (1/10 mm)	008	010	007
Cutting Length (mm)	4.2	4.2	2.5
Grit			
Shank Type	FG	FG	FG

Preparation Carbides



C-Series™					
Us No.	4C	35C	245C	330C	557C
Order No.	H1	H2	H245	H7	H31
Head Size (1/10 mm)	014	010	008	008	010
Cutting Length (mm)	—	0.9	2.8	2.0	3.9
Blade	—	—	—	—	—
Shank Type	FG	FG	FG	FG	FG

Finishing Diamonds



ET® Series								
Us No.	DET3F	DET6F	DET9F	—	—	DOS1F	—	—
Order No.	132F	134F	135F	8801	8888	8379	8368	8392
Head Size (1/10 mm)	008	014	014	023	012	023	016	016
Cutting Length (mm)	3.0	6.0	9.0	—	8.0	4.2	3.5	5.0
Grit								
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG

Finishing Carbides



ET® Series Safe-End										C-Series™	
Us No.	ET6	ET9	ET6F	ET6UF	ET9UF	OS1UF	7901	—	7408/OS1	7408CM	7901CM
Order No.	H134	H135	H134F	H134UF	H315UF	H379UF	H246	H274	H379	H379CM	H246CM
Head Size (1/10 mm)	014	014	014	014	014	023	009	016	023	023	009
Cutting Length (mm)	6.0	9.0	6.0	6.0	9.0	4.2	3.6	3.5	4.2	4.2	3.8
Blade											
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG



Finishing Strips

	Order No.	Thickness (mm)	Width (mm)	Length (mm)	Grit	Single/Double Sided
VisionFlex™	WS25F	0.10	2.5	147		Single
Gateway™	DS37EF	0.08	3.75	147		Single
Plastic Mylar	ETS1W	—	3.75	147	/	Single

Universal Composite Polishers



	BrioShine™ Single-Patient Sterile						
Order No.	SFLBC0-1	SBCPCOMP-1	SBPTCOMP-1	SBFLCP-C0	SFLBF-2	SBCPF-2	SBPTF-2
Head Size ($\frac{1}{10}$ mm)	100	060	040	055	100	060	040
Shape	Feather Lite™	Cup	Point	Feather Lite™ Cup	Feather Lite™	Cup	Point
Grit	Medium	Medium	Medium	Medium	Fine	Fine	Fine
Shank Type	RA	RA	RA	RA	RA	RA	RA



	Diacomp® Multi-Use			Composipro™	Winter Restorative Concepts		
Order No.	FLDC5M	FLDC5	DC3M.21	D209001	RWMLPT	RWFCUP	SPPP
Head Size ($\frac{1}{10}$ mm)	140	140	060	060	050	060	150
Cutting Length (mm)	1.6	1.6	7.5	5.5	10.0	7.5	2.5
Shape	Feather Lite™		Cup	Brush	Point	Cup	Knife Edge
Grit	Medium	Fine	Medium	Medium/Fine	Medium	Fine	Fine
Shank Type	RA	RA	RA	RA	RA	RA	RA



	EP™	Versaflex™
Order No.	EP1SM	VFD1LG
Head Size ($\frac{1}{10}$ mm)	80	120
Grit	Coarse	Super-Coarse
Type	Pop-on	Snap-on



Mandrels

Order No.	EPM21	310.21
Head Size ($\frac{1}{10}$ mm)	—	—
Type	Pop-on	Snap-on
Shank Type	RA	RA

INDIRECT RESTORATION

Brasseler's recommended diamonds, carbides, and polishers for each step of Crown & Bridge Restorations.

FULL CROWN

Porcelain Removal Diamonds



Order No.	5856	5837
Head Size (1/10 mm)	016	014
Cutting Length (mm)	8.0	8.0
Grit		
Shank Type	FG	FG

Zirconia & Lithium Disilicate Removal Diamonds



	DuraCut™		Dialite®	
Order No.	6856DC	TKS2	8856DF	
Head Size (1/10 mm)	016	014	016	
Cutting Length (mm)	8.0	8.0	8.0	
Grit				
Shank Type	FG	FG	FG	

Metal Cutting Carbides



Order No.	H34	H34L
Head Size (1/10 mm)	010	012
Cutting Length (mm)	2.0	3.5
Blade	—	—
Shank Type	RA	RA

Depth Cutting Diamonds



	Winter Restorative Concepts							
Order No.	WDLMIN	WDLCONV	WDLXT	WDLALC	WDL1.0	WDL1.5	WDL2.0	WDL2.5
Head Size (1/10 mm)	026	030	034	046	016	017	017	018
Depth Cut (mm)	0.3 / 0.5 / 0.7	0.5 / 0.7 / 0.9	0.8 / 1.0 / 1.2	1.0 / 1.2 / 1.5	1.0	1.5	2.0	2.5
Grit								
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG

Breaking Contact Diamonds



Order No.	KS0	5858	5850
Head Size (1/10 mm)	010	014	012
Cutting Length (mm)	7.0	8.0	10.0
Grit			
Shank Type	FG	FG	FG

Occlusal/Lingual Reduction Diamonds




Order No.	5379	5368	5811	909
Head Size (1/10 mm)	023	023	033	040
Cutting Length (mm)	4.2	5.0	4.2	1.0
Grit				
Shank Type	FG	FG	FG	FG

Preparation Diamonds



	Winter Restorative Concepts Multi-Radius											
Order No.	WS6M1.0	WS6F1.0	WS6M1.2	WS6F1.2	WS6M1.5	WS6F1.5	WS8M0.6	WS8F0.6	WS8M1.0	WS8F1.0	WS8M1.2	WS8F1.2
Head Size (1/10 mm)	016	016	018	018	020	020	014	014	018	018	020	020
Cutting Length (mm)	6.0	6.0	6.0	6.0	6.0	6.0	8.0	8.0	8.0	8.0	8.0	8.0
Tip Size (1/10 mm)	1.0	1.0	1.2	1.2	1.5	1.5	0.6	0.6	1.0	1.0	1.2	1.2
Grit												
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG	FG

Preparation Diamonds




KS Series									
Order No.	KS0	KS1	KS1L	KS2	KS3	KS4	KS5	KS6	KS7
Head Size ($\frac{1}{10}$ mm)	010	012	012	014	016	024	031	025	039
Cutting Length (mm)	7.0	7.0	9.0	8.0	8.0	4.2	6.2	9.0	1.0
Grit									
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG	FG


Preparation Diamonds



DuraCut™									
Order No.	878	878K	885	862	856	856	2856	6856DC	5856
Head Size ($\frac{1}{10}$ mm)	012	016	012	012	012	014	016	016	016
Cutting Length (mm)	8.0	8.0	8.0	8.0	8.0	8.0	8.0	8.0	8.0
Grit									
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG	FG




DuraBraz®									
Order No.	T5856	6856	BR6856	856	8856	856	6856	856	6856
Head Size ($\frac{1}{10}$ mm)	016	016	016	016	016	018	021	025	027
Cutting Length (mm)	8.0	8.0	8.0	8.0	8.0	8.0	8.0	8.0	8.0
Grit									
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG	FGSS



Order No.	6856	856L	856L	850	850	5847KR	8847KR	5847KR	8847KR
Head Size ($\frac{1}{10}$ mm)	029	016	018	016	018	016	016	018	018
Cutting Length (mm)	8.0	9.0	9.0	10.0	10.0	8.0	8.0	8.0	8.0
Grit									
Shank Type	FGSS	FG	FG	FG	FG	FG	FG	FG	FG

Margin Refinement Diamonds



Winter Restorative Concepts			
Order No.	WSMR0.6	WSMR1.0	WSMR1.4
Head Size ($\frac{1}{10}$ mm)	006	010	014
Cutting Length (mm)	0.17	0.27	0.25
Tip Size ($\frac{1}{10}$ mm)	0.6	1.0	1.4
Grit			
Shank Type	FG	FG	FG

INDIRECT RESTORATION



Brasseler's recommended diamonds, carbides, and polishers for each step of Crown, Bridge, and Veneer Restorations.

LAMINATE VENEER

Depth Cutting Diamonds

					
US No.	—	—	—	LVS1	LVS2
Order No.	828	828	828	834	834
Head Size ($\frac{1}{10}$ mm)	022	026	030	021	016
Cutting Length (mm)	1.2	1.2	1.2	6.0	6.0
Depth Cut (mm)	0.3	0.5	0.7	0.5	0.3
Grit					
Shank Type	FG	FG	FG	FG	FG

Preparation Diamonds

			
US No.	LVS3	LVS4	—
Order No.	6844	6844	5856
Head Size ($\frac{1}{10}$ mm)	016	014	016
Cutting Length (mm)	10.0/1.5	10.0/1.5	8.0
Depth Cut (mm)	—	—	—
Grit	Coarse/Fine	Coarse/Fine	
Shank Type	FG	FG	FG

Finishing Diamonds

				
	—	ET® Series		
US No.	—	DET6EF	DOS1F	DOS1EF
Order No.	8392	134EF	379F	379EF
Head Size ($\frac{1}{10}$ mm)	016	014	023	023
Cutting Length (mm)	5.0	6.0	4.2	4.2
Grit				
Shank Type	FG	FG	FG	FG

INLAY/ONLAY

Preparation Diamonds

				
Order No.	845KR	845KR	846KR	849L
Head Size ($\frac{1}{10}$ mm)	018	025	016	009
Cutting Length (mm)	4.0	4.0	6.0	5.0
Grit				
Shank Type	FG	FG	FG	FG

Finishing Diamond

	
Order No.	8274
Head Size ($\frac{1}{10}$ mm)	016
Cutting Length (mm)	3.5
Grit	
Shank Type	FG

Finishing Carbide

	
Order No.	H274UF
Head Size ($\frac{1}{10}$ mm)	016
Cutting Length (mm)	3.5
Blade	
Shank Type	FG

INTRA-ORAL ADJUSTING & POLISHING OF ALL CERAMIC MATERIALS

Adjusting Stone



Dialite® ZR

Order No.	LD11CZR
Head Size (1/10 mm)	030
Cutting Length (mm)	7.5
Shape	Point
Grit	Coarse
Shank Type	RA

Adjusting Diamonds



Dialite®

Order No.	8369DF	8856DF
Head Size (1/10 mm)	025	016
Cutting Length (mm)	5.5	8.0
Grit		
Shank Type	FG	FG

Zirconia Polishers



Dialite® ZR

Order No.	FLW14MZR	W17MZR	W16MZR	FLW14FZR	W17FZR	W16FZR
Head Size (1/10 mm)	140	060	040	140	060	040
Cutting Length (mm)	1.6	7.5	10.0	1.6	7.5	10.0
Shape	Feather Lite™	Cup	Point	Feather Lite™	Cup	Point
Grit	Medium	Medium	Medium	Fine	Fine	Fine
Shank Type	RA	RA	RA	RA	RA	RA

Lithium Disilicate Polishers



Dialite® LD

Order No.	FLW14MLD	W17MLD	W16MLD	FLW14FLD	W17FLD	W16FLD
Head Size (1/10 mm)	140	060	040	140	060	040
Cutting Length (mm)	1.6	7.5	10.0	1.6	7.5	10.0
Shape	Feather Lite™	Cup	Point	Feather Lite™	Cup	Point
Grit	Medium	Medium	Medium	Fine	Fine	Fine
Shank Type	RA	RA	RA	RA	RA	RA

Feldspathic Porcelain Polishers



Dialite® HP

Order No.	FLW14C	W17DCHP	FLW14M	W17DMHP	FLW14F	W17DFHP
Head Size (1/10 mm)	140	060	140	060	140	060
Cutting Length (mm)	1.6	7.5	1.6	7.5	1.6	7.5
Shape	Feather Lite™	Cup	Feather Lite™	Cup	Feather Lite™	Cup
Grit	Coarse	Coarse	Medium	Medium	Fine	Fine
Shank Type	RA	RA	RA	RA	RA	RA

Universal Ceramic Polishers



BrioShine™ Single-Patient Sterile

Order No.	SFLBCER-1	SBCPCER-1	SBPTCER-1	SFLBF-2	SBCPF-2	SBPTF-2
Head Size (1/10 mm)	100	060	040	100	060	040
Shape	Feather Lite™	Cup	Point	Feather Lite™	Cup	Point
Grit	Medium	Medium	Medium	Fine	Fine	Fine
Shank Type	RA	RA	RA	RA	RA	RA

INDIRECT RESTORATION

Brasseler's recommended diamonds, carbides, and polishers for each step of Crown, Bridge, and Veneer Restorations.


EXTRA-ORAL ADJUSTING & POLISHING OF ALL CERAMIC MAERIALS

Lithium Disilicate Grinder



Order No.	LD15M
Head Size (1/10 mm)	150
Cutting Length (mm)	3.0
Grit	Medium
Shape	Wheel
Shank Type	HP

Zirconia Polishers





Dialite® ZR				
Order No.	R17MZR	R17FZR	L26MZR	R20FZR
Head Size (1/10 mm)	170	170	260	200
Cutting Length (mm)	2.5	2.5	2.0	3.0
Grit	Medium	Fine	Medium	Fine
Shape	Wheel	Wheel	Knife Edge	Wheel
Shank Type	HP	HP	HP	HP

Lithium Disilicate Polishers



Dialite® LD		
Order No.	R17MLD	R17FLD
Head Size (1/10 mm)	170	170
Cutting Length (mm)	2.5	2.5
Grit	Medium	Fine
Shape	Wheel	Wheel
Shank Type	HP	HP


Feldspathic Porcelain Polishers



Dialite® HP				
Order No.	R17DCHP	R17DMHP	R17DFHP	FL26CHP
Head Size (1/10 mm)	170	170	170	260
Cutting Length (mm)	2.5	2.5	2.5	1.6
Grit	Coarse	Medium	Fine	Coarse
Shape	Wheel	Wheel	Wheel	Feather Lite™
Shank Type	HP	HP	HP	HP


PROVISIONALS

Diamond Discs



		HyperFlex™		VisionFlex™
Order No.	365	911HF	911H	934
Head Size (1/10 mm)	220	220	220	180
Coverage (mm)	3.0	3.0	3.0	3.0
Grit	Medium	Medium	Medium	Medium
Shank Type	HP	HP	HP	HP

Laboratory Carbide Cutters



	PSQ™	EF™		E™		SGE™	FSQ™
Order No.	H79PSQ	H257EF	H251EF	H78E	H351E	H79SGE	H79FSQ
Head Size (1/10 mm)	045	023	060	060	060	070	040
Cutting Length (mm)	13.0	13.0	14.0	12.0	11.0	14.0	13.0
Shank Type	HP	HP	HP	HP	HP	HP	HP



Polishers

	Scotch-Brite™	AcryliPro™	Winter Restorative Concepts	Dialite® PMMA	ET® ProviPro™
Order No.	25525HP	0674	RWGPP	PFL26FHP	SWR22F
Head Size (1/10 mm)	250	100	220	260	220
Cutting Length (mm)	—	24.0	3.0	1.6	3.0
Grit	Fine	Coarse	Medium	Fine	Fine
Shape	Brush	Bullet	Wheel	Feather Lite™	Wheel
Shank Type	HP	HP	Use with 305	HP	Use with 305



Mandrels


Order No.	305
Head Size (1/10 mm)	050
Shank Type	HP

SPECIALTY RESTORATIONS

Brasseler's recommended diamonds and carbides for Orthodontic, Surgical, and Endodontic procedures.

SPECIALTY

Orthodontic Debonding Carbides



Us No.	—	—	7675	7675RUF
Order No.	H48L	H48LUF	H375R	H375RU
Head Size (1/10 mm)	012	012	016	016
Cutting Length (mm)	8.0	8.0	8.0	8.0
Blade				
Shank Type	FG	FG	FG	FG

Surgical Carbides



Sabrecut®


Us No.	557	702	1558SC
Order No.	H31	H33	H031RS
Head Size (1/10 mm)	010	016	012
Cutting Length (mm)	4.2	4.4	4.0
Blade	—	—	—
Shank Type	FG SURG	FG SURG	FG SURG

Surgical Diamonds



Us No.	—
Order No.	859
Head Size (1/10 mm)	010
Cutting Length (mm)	11.0
Grit	
Shank Type	FG SURG

Bone Cutter Carbides



Us No.	—	—
Order No.	H161	H162
Head Size (1/10 mm)	016	016
Cutting Length (mm)	9.0	9.0
Blade	—	—
Shank Type	FG	FG

Endodontic Preparation Carbides



Sabrecut®

Us No.	1558SC
Order No.	H031RS
Head Size (1/10 mm)	012
Cutting Length (mm)	4.0
Blade	—
Shank Type	FG

Endodontic Preparation Diamonds



DuraCut™

Us No.	—
Order No.	6856DC
Head Size (1/10 mm)	016
Cutting Length (mm)	8.0
Grit	
Shank Type	FG

Endodontic Pilot Access Diamonds




DuraCut™

Dialite®

Us No.	—	—	—
Order No.	6801DC	15802	8801LDF
Head Size (1/10 mm)	014	018	014
Cutting Length (mm)	—	10.0	—
Grit			
Shank Type	FG	FG	FG


Endodontic Access Carbides



Sabrecut®

Us No.	—	1/2	4	6
Order No.	H269GK	H1	H1	H1
Head Size (1/10 mm)	016	006	014	018
Cutting Length (mm)	9.0	—	—	—
Shank Type	FG	FG SURG	FG SURG	FG SURG

Endodontic Access Diamonds




Us No.	—	—	—
Order No.	851	851	857
Head Size (1/10 mm)	012	016	014
Cutting Length (mm)	8.0	8.0	10.0
Grit			
Shank Type	FG	FG	FG

SINGLE-PATIENT DIAMONDS

Brasseler's recommended single-patient sterile diamonds for use in a variety of restorative procedures.

BrioPrep™ Diamonds



Us No.	—	—	—	—	—	BRIOLVS2
Order No.	BRI06830	BRI06835	BRI05858	BRI0859	BRI0801	BRI0834
Head Size (1/10 mm)	010	010	014	010	018	016
Cutting Length (mm)	2.8	4.0	8.0	11.0	—	6.0
Depth Cut (mm)	—	—	—	—	—	0.3
Grit						
Shank Type	FG	FG	FG	FG	FG	FG



Order No.	BRI05379	BRI08379	BRI08379	BRI05368	BRI08368	BRI05811	BRI05909
Head Size (1/10 mm)	023	018	023	023	023	037	040
Cutting Length (mm)	4.2	3.4	4.2	5.0	5.0	7.0	1.0
Grit							
Shank Type	FG	FG	FG	FG	FG	FG	FG




Us No.	—	—	—	—	—	BRIOKS0	BRIOKS1	BRIOKS2
Order No.	BRI05856	BRI05856	BRI05856	BRI05856L	BRI06850	BRI030005	BRI035005	BRI035005
Head Size (1/10 mm)	014	016	018	016	018	010	012	014
Cutting Length (mm)	8.0	8.0	8.0	9.0	10.0	7.0	7.0	8.0
Grit								
Shank Type	FG	FG	FG	FG	FG	FG	FG	FG




Order No.	BRI010839	BRI05847KR	BRI05847KR	BRI05862	BRI05878K
Head Size (1/10 mm)	014	016	018	012	016
Cutting Length (mm)	1.2	8.0	8.0	8.0	8.0
Grit					
Shank Type	FG	FG	FG	FG	FG


FRICTION GRIP



RIGHT-ANGLE LATCH



STRAIGHT HANDPIECE



TRIMMING & FINISHING CARBIDE BLADE CONFIGURATION

- 12 Blade ● 20 Blade ○ 30 Blade

ET® CARBIDE BLADE CONFIGURATION

- 8 Blade ● 15 Blade ● 16 Blade ○ 30 Blade


DIAMOND GRIT SIZES


- TM = 180 micron ● Medium = 100 micron
● DuraCut™ (DC) = 151 micron ● Fine = 30 micron
● Super-Coarse = 150 micron ● Extra-Fine = 15 micron
● Coarse = 125 micron ○ Ultra-Fine = 8 micron

HEAD/DIAMETER SIZES

Per ISO standards the diameter of a rotary instrument is defined as the widest part of the cutting surface.

1/10 mm	005	006	007	008	009	010	012	014	016	018	021	023	025	027	029
mm	0.5	0.6	0.7	0.8	0.9	1.0	1.2	1.4	1.6	1.8	2.1	2.3	2.5	2.7	2.9
inches	0.020	0.024	0.028	0.031	0.035	0.039	0.047	0.055	0.063	0.071	0.083	0.091	0.098	0.106	0.114
1/10 mm	031	033	035	037	040	042	045	047	050	055	060	065	070	075	080
mm	3.1	3.3	3.5	3.7	4.0	4.2	4.5	4.7	5.0	5.5	6.0	6.5	7.0	7.5	8.0
inches	0.122	0.130	0.138	0.148	0.157	0.165	0.177	0.185	0.197	0.217	0.236	0.256	0.276	0.295	0.315






DIAMONDS



CARBIDES



POLISHERS



PROCEDURE SYSTEMS



BIOCERAMICS



Since 1976, our absolute focus has been to develop products and provide services that support the practice of core operator procedures. Today, Brasseler® is the premier dental instrumentation company in North America offering the most extensive selection of dental and surgical instrumentation under one brand.



HANDPIECES



ENDODONTICS



SMALL EQUIPMENT



HYGIENE



HAND INSTRUMENTS



ORAL SURGERY



LABORATORY INSTRUMENTS



Visit our website at BrasselerCanada.com
To order call 800.363.3838

