



# EXCLUSIVE GPO STARTER PACKAGE PROMOTIONS



## Torq Control® Universal Torque Wrench

PART NUMBER: 15501 | LIST PRICE: \$799

Buy 1 Torq Control Wrench at 20% discount | **\$639**  
Instaquote Bundle Code: **NAM-ANTH-005**

Buy 2 Torq Control Wrenches at 35% discount | **\$1,039**  
Instaquote Bundle Code: **NAM-ANTH-06**


Torq Control® is a universal torque wrench that is designed to ensure accurate tightening regardless of the type of implant connection or accessibility challenges.

Tightening accuracy is crucial for successful implant treatment and Torq Control® has been specifically designed to meet these needs for all prosthetic procedures, with added control through automatic declutching.

The original Torq Control is a must-have, especially for full-arch restorations.

**ADAPTED FOR  
ALL PROTOCOLS**

7 torque values



**100° MICRO-HEAD**

Posterior sector easy access  
+ Excellent visibility

**UNIVERSAL  
CONNECTION**

All prosthetic mandrels  
Any implant system

**ERGONOMIC DESIGN**


Patient comfort  
Ease of use - 130 g only

**MANUAL INSTRUMENT**

Screw-in feel for full control

**AUTOMATIC  
DECLUTCHING**

Precision + Control



**Synergy**  
dental partners

## COMPATIBLE WITH ALL IMPLANT SYSTEMS!

- Universal connection for all dental implant systems
- 7 torque values, from 10 to 35 Ncm, adapted for all protocols

## DESIGNED FOR EFFICIENCY

- Micro-head, 100° angulated for easy access and optimal visibility, even in posterior sectors

## MAXIMUM COMFORT OF USE

- Manipulation 100% outside of mouth
- No lateral movement

## SO EASY TO HANDLE

- Only 130 g for a better freedom of movement

## CONTROL FEATURES

- High precision torque tightening
- Automatic declutching when the selected torque is reached

## EASY MAINTENANCE

- One-piece design
- Sterilizable in autoclave, 135°

## GUARANTEE OF QUALITY

- Made in the Mont-Blanc Valley, France, since 1997
- Top precision with more than 50 stages of high-quality manufacturing



### CALIBRATION

In order to maintain the precision of the torque values, it is recommended that the torque wrench be returned for calibration every 2,000 uses or every 2 years.

Suggested Third Party Calibration  
**Xhandpiece Systems** / (800) 319-3243 /  
xhandpiecesystems.com

**IN-CAL** / (603) 883-5558 / in-cal.com

**CANADA / Marsh Metrology**  
1 (905) 331-9783 / marshmetrology.com



### Straumann North American Headquarters

Straumann USA, LLC  
60 Minuteman Road  
Andover, MA 01810  
Phone 800/448 8168 (US) • 800/363 4024 (CA)  
Fax 978/747 2490  
www.straumann.us • www.straumann.ca

ifu.anthogyr.com

© Straumann USA, LLC 2025. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.