



MULTI-PURPOSE. TEMPORARY BITE PLANE

DESIGNED FOR SHORT-TERM OVERNIGHT WEAR

Easy to customize in minutes with Futar or Panasil, QuickSplint is a practical and affordable anterior bite plane. Auxiliary staff can fabricate QuickSplint for immediate delivery and you can reline the QuickSplint for staged restorative procedures.



Fits upper or lower arch.



QuickSplint®keeps the posterior teeth out of occlusion. This inhibits cleching and reduces bite force, which is beneficial in many applications:

PROTECT TEMPORARIES AND NEWLY PLACED VENEERS

Use QuickSplint®to inhibit overnight clenching, ensuring that patients will not dislodge temporaries or veneers. Also helps patients habituate to occlusal changes.

SIMPLIFY DEPROGRAMMING FOR CENTRIC BITE RECORDS

Use QuickSplint to inhibit overnight clenching, ensuring that patients will not dislodge temporaries or veneers. Also helps patients habituate to occlusal changes.

SIMPLIFY DEPROGRAMMING FOR CENTRIC BITE RECORDS

QuickSplint®levels the occlusal plane and seats the condyles, allowing you to take repeatable bite records.



FILL



ALIGN



TRIM









DIAGNOSE SEVERITY OF BRUXISM

Provide QuickSplint on a short-term basis to communicate with patients and raise awareness of the seriousness of their bruxism.

This is ideal for case planning and informed consents for cosmetic and restorative patients.

ASSIST WITH CLEAR ALIGNER TRAYS

Use QuickSplint over the trays to interfere with clenching behavior at night. Reline the QuickSplint for changes in aligners.

PROMOTE POST-OP HEALING

Deliver QuickSplint chairside on upper or lower arch, to keep teeth out of occlusion during the healing phase of perio, endo, extractions and implant procedures to minimize microtrauma.



DELIVER

YOUR QUICKSPLINT® PARTNERSHIP PRICING

	Retail Price	Member Price
QuickSplint®(12) Pack	\$240 ^{.00}	\$192 ^{.00}
QuickSplint®(100) Pack	\$1,500 .00	\$1,200 ^{.00}



Call Kettenbach Direct! (877) 532-2123 www.kettenbachusa.com