

Direct Restorative

Venus® Pure

Simplify your direct restorative workflow.

Venus Diamond PURE | Venus Pearl PURE

Four shades to seamlessly blend with surrounding tooth structure, Venus PURE resins optimize workflow efficiency without compromising on strength and durability.

- Easy to use; a perfect fit for 90% of daily restorative procedures
- Monochromatic layering to bring simplicity to everyday anterior and posterior restorations

Venus Pearl PURE Shades – Refills	Darby Item #
Venus Pearl PLT Refill 20x0,2G Light	525-6350
Venus Pearl PLT Refill 20x0,2G Medium	525-6351
Venus Pearl PLT Refill 20x0,2G Dark	525-6352
Venus Pearl PLT Refill 20x0,2G Bleach	525-6353
Venus Pearl Syringe Refill 1x3G Light	525-6354
Venus Pearl Syringe Refill 1x3G Medium	525-6355
Venus Pearl Syringe Refill 1x3G Dark	525-6356
Venus Pearl Syringe Refill 1x3G Bleach	525-6357

Venus Diamond PURE Shades – Refills	Darby Item #
Venus Diamond PLT Refill 20x0,25G Light	525-6340
Venus Diamond PLT Refill 20x0,25G Medium	525-6341
Venus Diamond PLT Refill 20x0,25G Dark	525-6342
Venus Diamond PLT Refill 20x0,25G Bleach	525-6343
Venus Diamond Syringe Refill 1x4G Light	525-6344
Venus Diamond Syringe Refill 1x4G Medium	525-6345
Venus Diamond Syringe Refill 1x4G Dark	525-6346
Venus Diamond Syringe Refill 1x4G Bleach	525-6347

NEW!

Synergy
dental partners

darby



KULZER
MITSUI CHEMICALS GROUP



Venus Pearl PURE Shades – Kits	Darby Item #
Venus Pearl Pure Kit PLT 20x0.2g PLTs each of Venus Pearl shades Light + Medium; 2g syringe of Venus Bulk Flow ONE; 2ml bottle of iBOND Universal; 5 Flow Cannulas; Translation card; Pictorial card iBOND Universal; IFUs	525-6359
Venus Pearl Pure Kit Syringe 3g syringes each of Pearl shades Light + Medium; 2g syringe of Venus Bulk Flow ONE; 2ml bottle of iBOND Universal; 5 Flow Cannulas; Translation card; Pictorial card iBOND Universal; IFUs	525-6358

Venus Diamond PURE Shades – Kits	Darby Item #
Venus Diamond Pure Kit PLT 20x0.25g PLTs each of Venus Diamond shades Light + Medium; 2g syringe of Venus Bulk Flow ONE; 2ml bottle of iBOND Universal; 5 Flow Cannulas; Translation card; Pictorial card iBOND Universal; IFUs	525-6348
Venus Diamond Pure Kit Syringe 4g syringes each of Venus Diamond shades Light + Medium; 2g syringe of Venus Bulk Flow ONE; 2ml bottle of iBOND Universal; 5 Flow Cannulas; Translation card; Pictorial card iBOND Universal; IFUs	525-6349